

# PE-04ATA50-P

## DYRK1A (316-330) pY319+pY321+pS324+pY327 Peptide Powder

15-mer immunogen and phosphatase substrate phosphopeptide based on DYRK1A



# KINEXUS

Address: 8755 Ash Street, Suite 1  
Vancouver, British Columbia,  
Canada V6P 6T3

Email: [info@kinexus.ca](mailto:info@kinexus.ca)  
Phone: 604-323-2547

### Target Protein

<b>Name Long:</b>	Dual specificity tyrosine-phosphorylation-regulated kinase 1A
<b>Name Alias:</b>	DYRK1A; DYRK; MNB; MNBH; HP86; DYRK1; CCDS13653.1; ENSG00000157540
<b>Species Origin:</b>	Human
<b>UniProt ID:</b>	Q13627

### Peptide Structure

<b>Peptide Name:</b>	DYRK1A (316-330) pY319+pY321+pS324+pY327
<b>Peptide Origin:</b>	In the protein kinase catalytic domain activation T loop region between subdomains VII and VIII.
<b>Peptide Sequence Location:</b>	Q316-P330
<b>Peptide Sequence:</b>	QRI(pY)Q(pY)IQ(pS)RF(pY)RSP
<b>Peptide N-Terminus:</b>	Free amino
<b>Peptide C-Terminus:</b>	Amide
<b>Peptide Modifications Other:</b>	Phosphorylated

### Production

<b>Peptide Production Method:</b>	Solid-phase peptide synthesis
<b>Calculated Peptide Mass:</b>	2498.6
<b>Observed Peptide Mass:</b>	2495.8
<b>% Peptide Purity:</b>	50.7
<b>Peptide Appearance:</b>	White powder
<b>Peptide Form:</b>	Solid
<b>Peptide Solubility:</b>	Dissolve in 50 µl DMSO and dilute to desired concentration with water or aqueous buffer
<b>Lot Number:</b>	KLP04CAB-07
<b>Amount:</b>	1 mg
<b>Storage Conditions:</b>	Frozen at -20 °C
<b>Storage Stability:</b>	Over 1 year at -20 °C

### Applications

<b>Product Use:</b>	This phosphopeptide may be useful as a substrate for screening the phosphatase activity of protein phosphatases.
---------------------	--

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit [www.kinexusproducts.ca](http://www.kinexusproducts.ca) or contact us at 1-866-KINASES (546-2737)